



8 in



# SM-500L

8 inch color high resolution TFT display

Portability, suitable for at patient's bedside and transferring monitoring

Unique blood pressure measurement technique

Anti high-frequency electro-tome and defibrillation, digital filter

Built-in rechargeable battery, DC/AC



# SM-500L

8inchTFT



## Picture:

- Big font display interface
- Trend graph display interface
- Standard display interface
- ECG, NIBP, SPO2, TEMP interface
- ECG full lead display interface
- Built-in rechargeable battery and recorder

## Feature:

- Display with 9 waveform
- 720 hours trend data
- 1000 NIBP data storage
- 200 alarm event data storage
- Sound, light, human voice alarm
- 12 hours waveform review
- Wire and wireless network
- Nurse call
- analog output
- Support software upgrade online
- Arrhythmia analysis, ST segment analysis
- Support touch screen

## Standard:

ECG, RESP, NIBP, SpO2, TEMP

## Optional:

Recorder, IBP,touch screen



## Sonostar Technologies Co., Limited

Address: #31, Shuiwan road, Jida, Zhuhai, Guangdong, China  
Tel: +86 13168355297 MSN: sonostar@163.com Skype: sonostar-ut  
Web: www.sonostar.net Mail: sonostar@sonostar.net

